

Model 7407 Vibrating Sample Magnetometer

Introduction

Lake Shore's Vibrating Sample Magnetometers perform magnetic measurements for materials research and development, quality control, and production testing.

The Model 7407 is capable of characterizing a variety of particulate and continuous magnetic media materials including; audio, video, and digital data tapes, flexible media, magneto-optical materials, sputtered and plated thin film materials including multilayer GMR, CMR, exchange-bias and spin valve materials.



Measurements

The following parameters are either measured directly or can easily be derived through the software

- Hysteresis Loops
 - Saturation magnetization (M_{SAT}), Retentivity or remanent magnetization (M_{REM})
 - Coercivity (H_C), S^* , Slope at H_C , value of dM/dH or differential susceptibility at H_C
 - Switching field distribution (SFD)
 - Flatness, Squareness ratio (SQR)
- Minor hysteresis loops
- Initial magnetization curve
- DC remanence
- AC remanence
- Vector measurements (m_x and m_y)
- Magnetization data as a function of time

Materials

All types of magnetic materials:

- Diamagnetic, Paramagnetic, Ferromagnetic, Ferrimagnetic, Antiferromagnetic materials and Anisotropic materials
- Particulate and continuous magnetic recording materials and GMR, CMR, exchange biased and spin-valve materials
- Magnetic-optical materials
- Bulk materials, powders, thin films, single crystals, and liquids are readily accommodated

Features

- Computer-automated data collection system providing field strengths to 23 kG
- Measurement of moments as small as 0.1×10^{-6} emu in room temperature magnetic fields ranging from -23 to +23 kG
- Adjustable magnet air gap permits adjusting magnet/coils to suit sample and field strength requirements
- Bipolar power supply provides smooth continuous transition through zero field
- Fast data acquisition - average sample run (hysteresis loop) over full field range typically requires only minutes
- Windows™ NT/2000 menu driven color graphic software for system operation, data acquisition, and analysis. System software includes operation and control of the magnet power supply, VSM control unit, and gaussmeter. Real-time feedback of processed magnetic moment measurement data can be displayed in either graphical or tabular format.
- Water cooled magnet coils provide excellent field stability when high power is required to achieve the maximum field capability

System Specifications

General

Magnetic Moment

| | |
|--------------------------------|--|
| Dynamic range | 0.1 x 10 ⁻⁶ emu to 10 ³ emu (room temperature) |
| Time constants | 0.1 second, 0.3 second, 1.0 second, 3.0 seconds, or 10.0 seconds |
| Output stability | Better than ±0.05% of full scale per day for fixed coil geometry at constant field and temperature |
| Absolute accuracy | Better than 2% of reading ±0.2% of full scale |
| Reproducibility | Better than ±1%, or ±0.15% of full scale, whichever is greater |
| Field accuracy in gauss | 1% of reading or ±0.05% of full scale |

Model 7407 Performance Specifications

| | | | | |
|--------------------------------|--|-----------------------------|-----------------------------|-----------------------------|
| Magnet pole face | 7.6 cm (3") | 15.2 cm (6") | | |
| Coil set | 730EMSC - 1" coils | 730ESC - 2" coils | | |
| Operating air gap range | Min 1.6 cm (0.63") | Max 4.1 cm (1.6") | Min 3.8 cm (1.5") | Max 5.1 cm (2.0") |
| Sample access | 0.33 cm (0.13") | 2.5 cm (1.0") | 1.3 cm (0.5") | 2.5 cm (1.0") |
| Magnetic field | 23 kG | 15 kG | 12 kG | 10 kG |
| Noise/sensitivity | 0.1 μemu at 0.63" air gap and room temperature | | | |

Equipment

Lake Shore Model 7407 consisting of:

Model 736 VSM Electronics with integrated gaussmeter

Resolution ± 1 part out of 300,000

Precision Up to 0.0007% of full scale for 350 G and above ranges

Model 730ESC pick-up coils or Model 730EMSC mini pick-up coils (Select one only)

VSM head assembly and mounting structure

Instrumentation console

Electromagnet Model EM7-HV Variable Gap Electromagnet Includes one set of pole caps with 3" or 6" pole diameter

Magnet pole face diameter 7.62 cm (3") or 15.24 cm (6")

Pole diameter 17.8 cm (7")

Flow rate 3 gallons per minute (11 liters per min.)

Cooling water requirements Tap water or closed cooling system (optional chiller available)

Model 665 Bipolar Power Supply

Maximum Output ± 50 V/ ± 100 A (5 kW)

AC line input 3-phase and neutral 208, 220, 380, 400 Vrms at 50/60 Hz

Flow rate: 2.2 gallons per minute (8 liters per minute)

Cooling water requirements: Tap water or closed cooling system (optional chiller available)

+11 °C to +25 °C

Sample holders (3 of each holder supplied)

Powder 3 mm (0.13") diameter x 6.35 mm (0.25") long

Thin-film bottom 6.35 mm (0.25") diameter

Thin-film side 6.35 mm (0.25") long x 12.7 mm (0.5") wide

Liquid (holder purchased separately) 3 mm (0.13") diameter x 6.35 mm (0.25") long; nominal volume 100 μ l

Computer

700 MHz Intel processor, 64 Meg of RAM, hard drive, 3.5" 1.44 MB floppy drive, CD-ROM, 17" SVGA color monitor with 1024x768 resolution, Windows™ NT/2000, and National Instruments PCI-GPIB IEEE-488 interface.

IDEAS™ VSM Software

Windows™ NT/2000 menu driven, enhanced color-graphic software for system operation, data acquisition and analysis. System software includes operation and control of the magnet power supply, VSM control unit, and gaussmeter. Real-time feedback of processed magnetic moment measurement data can be displayed in either graphical or tabular format.

Printer HP InkJet printer

Shipping weight Four (4) crates totaling 1093 kg (2405 lbs)

Copyright Lake Shore 07/02